



PuroSPIN™ Fibrous Tissue RNA Purification Kit

- *RNA Purification from Fibrous Animal Tissue*
- *RNA Cleanup*

#NK053-50, #NK053-250

Product Manual

Version 2.1

www.lunanano.com

Kit Contents

Component	50 Preps (#NK053-50)	250 Preps (#NK053-250)
Lysis Buffer LB-R	40 mL	200 mL
Wash Buffer WB-R1	50 mL	225 mL
Wash Buffer WB-R2	60 mL *	225 mL *
Proteinase K (20 mg/mL), Nuclease-free	600 μ L	3 x 1 mL
DNase I, RNase-free	550 U	5 x 550 U
DNase I Reaction Buffer	5 mL	25 mL
Elution Buffer EB	15 mL	50 mL
Water, Nuclease-free	15 mL	50 mL
PuroSPIN™ MINI Spin Columns	50	250
Collection Tubes (2 mL)	50	250

* After addition of 96-100% Ethanol

Storage

DNase I and **Proteinase K** need to be stored at $\leq -20^{\circ}\text{C}$. Once reconstituted in water, store **DNase I** at 4°C (< 6 weeks), or in one-time-use aliquots at $\leq -20^{\circ}\text{C}$. All other components can be stored at room temperature for at least 1 year.

Additional materials and equipment required

- 2 M DTT or 14 M β -mercaptoethanol
- 96-100% Ethanol
- Pipettes, sterile nuclease-free tips
- 1.5 mL RNase-free microcentrifuge tubes
- Microcentrifuge (with 12,000g capability)
- Vortex (if not available, can be replaced with vigorous pipetting)
- Appropriate personal protective equipment (gloves, lab coat)
- Water bath or heating block (heated to 55°C).
- Rotor-stator homogenizer, or bead-beater apparatus (with 2-5 mm zirconia or steel beads), or mortar and pestle with liquid nitrogen for tissue digestion


Kit Description

PuroSPIN™ Fibrous Tissue RNA Purification Kit is designed for extraction of total RNA from fibrous animal tissue, such as skin, skeletal muscle, lung, or heart. It greatly simplifies the extraction process by incorporating a convenient easy-to-use extraction technology based on nucleic acid binding silica spin-columns. The tissue is first disrupted through homogenization and Proteinase K treatment. The RNA is then released from the cells by the treatment with the **Lysis Buffer LB-R**. Following the lysis step, RNA is loaded onto silica-based spin columns and washed to remove protein, salt and other contaminants. ***The kit includes on-column DNase I digestion step to remove any contaminating DNA.*** Finally, purified RNA is eluted from the column into a clean microcentrifuge tube in **Elution Buffer EB**, or **nuclease-free water**.

The kit can also be used to perform RNA cleanup after DNase I or other enzymatic treatments. High quality RNA extracted with the **PuroSPIN™ Fibrous Tissue RNA Purification Kit** can be used for a range of downstream applications, such as reverse-transcription PCR or qPCR, northern blotting, microarray analysis.

For purification of RNA from soft tissue, such as spleen, kidney, liver, or brain, or from bacteria, yeast, or mammalian cells, please use our **PuroSPIN™ Total RNA Purification Kit** (#NK051-50, #NK051-250).

Experimental Workflow



Important Notes

BEFORE USE:

- Add 96-100% Ethanol to the **Wash Buffer WB-R2**:
 - **#NK053-50** (50 preps): add **45 mL** of Ethanol
 - **#NK053-250** (250 preps): add **170 mL** of Ethanol
- Reconstitute each tube of lyophilized **DNase I** in **550 μ L** of nuclease-free water
 - Pulse-spin the tube before water addition to ensure the pellet is at the bottom.
 - **Do not vortex**. Vortexing might cause DNase I denaturation. Flip the tube upside down a few times to reconstitute the enzyme.
 - Store **DNase I** at **4°C** (< 6 weeks), or in one-time-use aliquots at \leq **-20°C**.
 - **#NK053-250**: reconstitute one tube at a time to prevent **DNase I** degradation.



Lysis Buffer LB-R contains **guanidine hydrochloride**. It is a known irritant that is harmful if inhaled, swallowed, or contacted with skin. Wear appropriate protective equipment when working with this reagent.

Within 30 minutes of starting RNA purification, make fresh **Working Buffer LB-R** by adding 20 μ L of 2 M DTT or 14 M β -mercaptoethanol solution to each 1 mL of **Lysis Buffer LB-R** used.

If precipitate is observed in the **Lysis Buffer LB-R** or **Wash Buffer WB-R1**, redissolve it by warming the solution to 37°C. Cool the solution back to room temperature before use.

RNA is highly prone to degradation by the RNase enzymes. RNases are very stable enzymes that are present on human skin and commonly contaminate lab equipment and surfaces. To minimize RNase contamination always wear gloves when performing the extraction, use RNase-free filtered tips, and pre-treat all non-disposable equipment (pipettes, vortex, centrifuge) and surfaces with RNase removal solution. Ensure all of the buffers are tightly closed when not in use. Take particular care with the **Elution Buffer EB** and **nuclease-free water** solutions, since these are used for storage of eluted RNA.

In addition to possibility of degradation by contaminating RNase enzymes, RNA can also spontaneously break down through auto-hydrolysis. To minimize the probability of both hydrolysis and RNase degradation, the samples should either be kept on ice or frozen prior to RNA extraction. Both bacterial and mammalian cells can be pelleted and stored at -70°C until the time of RNA extraction. Lysis and homogenization steps should be performed as soon as possible after the sample is removed from ice or thawed. If fresh samples are used, the time between sample collection and RNA extraction should be minimized. Following the extraction, RNA should either be (i) kept on ice and used as soon as possible for the downstream application, or (ii) immediately frozen and stored at -70°C until needed.

Protocol A. Purification of Total RNA from Fibrous Animal Tissue

Within 30 min of starting the RNA extraction, prepare the following reagents:

- Make **Working Buffer LB-R** by adding 20 μL of 2 M DTT or 14 M β -mercaptoethanol solution to each 1 mL of **Lysis Buffer LB-R** used. Each 1 mL of **Working Buffer LB-R** is sufficient for 2 reactions.

Step	Procedure
1	Collect up to 30 mg of fresh or frozen tissue.
2	<p>Homogenize the tissue using one of the following methods:</p> <p>a) <u>Rotor-stator homogenizer</u></p> <ul style="list-style-type: none">▪ Transfer the tissue sample into an appropriate container. Add 500 μL of Working Buffer LB-R into the tube.▪ Homogenize the sample for 20-50 seconds using standard rotor-stator homogenizer (please refer to the apparatus manufacturer's manual for specific instructions). <p>b) <u>Bead-beater homogenizer</u></p> <ul style="list-style-type: none">▪ Cut the tissue into small pieces ~ 1 mm in cross-section. Transfer the samples into a new RNase-free 2 mL screw cap tube (ensure the tube is compatible with the bead-beater apparatus).▪ Add 500 μL of Working Buffer LB-R into the tube.▪ Add the beads to the tube<ul style="list-style-type: none">○ For fibrous tissue such as muscle/heart/lung, use ~ 10 of 2.8 mm steel beads○ For elastic tissue such as skin/sclera/cartilage use 1 of 6.4 mm steel bead▪ Disrupt the tissue for 2-3 min using a bead-beater homogenizer (please refer to bead-beater apparatus manufacturer's manual for specific instructions).▪ Centrifuge the samples at 100g for 1 min to pellet the beads. Transfer the supernatant into a new RNase-free 1.5 mL microcentrifuge tube. <p>c) <u>Mortar and pestle</u></p> <p>Freeze the tissue by immersing the tube containing it into liquid nitrogen. Grind the tissue thoroughly using mortar and pestle. Add 500 μL of Working Buffer LB-R to the tissue and disrupt it further by passing it 5-6 times through a 20-gauge needle using RNase-free syringe.</p>
3	Add 300 μL of nuclease-free water to the homogenized sample.
4	Add 10 μL of Proteinase K solution to the sample. Mix well by pipetting.

5	Incubate the samples for 10 min at 55°C.
6	Centrifuge the tube at 12,000g for 3 min. Carefully collect the supernatant and transfer it into a new RNase-free 1.5 mL microcentrifuge tube.
7	Add 500 µL of 96-100% Ethanol. Invert a few times to mix.
8	Transfer 700 µL of the solution onto the PuroSPIN™ MINI Spin Column inserted into the provided 2 mL collection tube. Centrifuge the column at 12,000g for 2 min to load the nucleic acids onto the column. Discard the flow-through.
9	Repeat <u>Step 8</u> until all of the sample has been loaded onto the column.
10	Add 350 µL of Wash Buffer WB-R1 to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
11	Centrifuge at 12,000g for 1 min to pre-wash and dry out the membrane.
12	Mix 10 µL of DNase I with 70 µL of DNase I Reaction Buffer .
13	Add 80 µL of buffered DNase I mixture solution to the column.
14	Incubate the column for 15 min at room temperature to allow DNA digestion to occur.
15	Add 350 µL of Wash Buffer WB-R1 to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
16	Add 700 µL of Wash Buffer WB-R2 (<i>with Ethanol added</i>) to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
17	Centrifuge at 12,000g for 1 min to dry out the membrane.
18	Transfer the PuroSPIN™ MINI Spin Column into a new RNase-free 1.5 mL microcentrifuge tube.
19	<p>Pipette 50-100 µL of Elution Buffer EB or nuclease-free water onto the center of the column. Centrifuge at 12,000g for 1 min to elute the RNA. Keep the flow-through that contains the product and discard the column.</p> <p><i>NOTE: Use lower elution volume to maximize the concentration of extracted RNA and larger elution volume to increase the total RNA recovery. To maximize total RNA recovery, repeat the elution with additional 100 µL of included Elution Buffer EB or nuclease-free water.</i></p>

Protocol B. RNA Cleanup Protocol

PuroSPIN™ Fibrous Tissue RNA Purification Kit can be used for cleanup of RNA preparation following enzymatic reactions, such as DNase I treatment (**Protocol D** of this kit). If RNA solution contains contaminating DNA, please follow **Protocol C**.

NOTE: For this protocol use **Lysis Buffer LB-R** directly without the addition of DTT or β -mercaptoethanol.

Step	Procedure
1	Adjust the volume of the RNA mixture to ~ 100 μ L with nuclease-free water .
2	Add 250 μ L of the Lysis Buffer LB-R (<i>no β-mercaptoethanol or DTT added</i>). Mix by vortexing or pipetting.
3	Add 350 μ L of 96-100% Ethanol. Invert a few times to mix.
4	Transfer 700 μ L of the solution onto the PuroSPIN™ MINI Spin Column inserted into the provided 2 mL collection tube. Centrifuge the column at 12,000g for 2 min to load the nucleic acids onto the column. Discard the flow through.
5	Add 700 μ L of Wash Buffer WB-R1 to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
6	Add 700 μ L of Wash Buffer WB-R2 (<i>with Ethanol added</i>) to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
7	Centrifuge at 12,000g for 1 min to dry out the membrane.
8	Transfer the PuroSPIN™ MINI Spin Column into a new RNase-free 1.5 mL microcentrifuge tube.
9	Pipette 50-100 μ L of included Elution Buffer EB or nuclease-free water onto the center of the column. Centrifuge at 12,000g for 1 min to elute the RNA. Keep the flow-through that contains the product and discard the column. <i>NOTE: Use lower elution volume to maximize the concentration of extracted RNA and larger elution volume to maximize the total RNA recovery. To maximize total RNA recovery, repeat the elution with additional 100 μL of included Elution Buffer EB or nuclease-free water.</i>

Protocol C. RNA Cleanup with DNA Removal Protocol

This protocol includes a DNase I on-column digestion step designed to remove any DNA contaminating the RNA solution.

NOTE: For this protocol use **Lysis Buffer LB-R** directly without the addition of DTT or β -mercaptoethanol.

Step	Procedure
1	Adjust the volume of the RNA mixture to ~ 100 μ L with nuclease-free water .
2	Add 250 μ L of the Lysis Buffer LB-R (<i>no or β-mercaptoethanol or DTT added</i>). Mix by vortexing or pipetting.
3	Add 350 μ L of 96-100% Ethanol. Invert a few times to mix.
4	Transfer 700 μ L of the solution into the PuroSPIN™ MINI Spin Column inserted into the 2 mL collection tube. Centrifuge the column at 12,000g for 2 min to load the nucleic acids onto the column. Discard the flow through.
5	Add 350 μ L of Wash Buffer WB-R1 to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
6	Centrifuge at 12,000g for 1 min to dry out the membrane.
7	Mix 10 μ L of DNase I with 70 μ L of DNase I Reaction Buffer .
8	Add 80 μ L of the DNase I mixture solution to the column.
9	Incubate the column for 15 min at room temperature to allow DNA digestion to occur.
10	Add 350 μ L of Wash Buffer WB-R1 to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
11	Add 700 μ L of Wash Buffer WB-R2 (<i>with Ethanol added</i>) to the column. Centrifuge at 12,000g for 1 min. Discard the flow-through.
12	Centrifuge at 12,000g for 1 min to dry out the membrane.
13	Transfer the PuroSPIN™ MINI Spin Column into a new RNase-free 1.5 mL microcentrifuge tube.

14	<p>Pipette 50-100 μL of included Elution Buffer EB or nuclease-free water onto the center of the column. Centrifuge at 12,000g for 1 min to elute the RNA. Keep the flow-through that contains the product and discard the column.</p> <p><i>NOTE: Use lower elution volume to maximize the concentration of extracted RNA and larger elution volume to maximize the total RNA recovery. To maximize total RNA recovery, repeat the elution with additional 100 μL of included Elution Buffer EB or nuclease-free water.</i></p>
----	--

Protocol D. Genomic DNA Removal from RNA Preparation

The protocol describes the in-solution removal of DNA from RNA preparation by **DNase I** treatment. The protocol is provided for 1 μ g of RNA. It can be scaled up to the appropriate amount of RNA used.

NOTE: 10X DNase Reaction Buffer (100 mM Tris-HCl, pH 7.5, 25 mM MgCl₂, 1 mM CaCl₂) is not provided with the kit.


Step	Procedure
1	<p>Set up the following reaction:</p> <ul style="list-style-type: none"> • RNA Preparation..... 1 μg • DNase I, RNase-free..... 1 unit • 10X DNase Reaction Buffer... 1 μL • RNase free water..... to 10 μL
2	Incubate at 37°C for 30 min.
3	<p>The reaction can be terminated either by thermal inactivation of DNase I, or by performing RNA purification.</p> <p>a) <u>Inactivation</u> Add 1 μL of 50 mM EDTA. Incubate at 65°C for 10 min. EDTA is added to chelate divalent cations to prevent RNA hydrolysis during the heating step.</p> <p>b) <u>Purification</u> Purify DNase I treated RNA using</p> <p><u>Protocol B</u> of this kit.</p>

Troubleshooting

Problem	Possible causes and solutions
Column clogging	<p>Too much starting material was used. Reduce the amount of material used for the lysate preparation. Do not use higher amount of material than indicated in the protocol.</p> <p>Starting material was not completely disrupted. Perform more thorough homogenization. Use smaller amount of starting material.</p> <p>Cellular debris was not completely removed. Perform centrifugation to remove cellular debris as indicated in the protocol.</p>
Low RNA yield	<p>Ethanol was not added to Wash Buffer WB2 Add the indicated amount of 96-100%.</p> <p>Too much starting material was used. Reduce the amount of material used for the lysate preparation. Do not use higher amount of material than indicated in the protocol.</p> <p>Starting material was not completely disrupted. Perform more thorough homogenization. Use smaller amount of starting material.</p> <p>Ethanol was not added to the lysate. Make sure the indicated amount of ethanol is added to the sample lysate before loading it onto the column.</p>
Degraded RNA	<p>RNase contamination Use RNase-free tips and tubes. Clean work surfaces and non-disposable items with RNase removal solution. Wear gloves.</p> <p>Inappropriate handling of the sample. When RNA is extracted from fresh samples, the samples should be kept on ice, lysis and homogenization should be performed as soon as possible. Frozen samples should be thawed immediately prior to the lysis step. Once extracted, RNA should be immediately placed on ice and stored between -20°C and -70°C.</p>
DNA contamination	<p>On-column DNase digestion was not performed. Follow the protocol steps to perform on-column DNase digestion.</p>

Notes

Safety

Reagent	Risk Symbols	Risk and Safety Phrases
Lysis Buffer LB-R		<p>Contains <i>guanidine hydrochloride</i></p> <ul style="list-style-type: none"> • Harmful if swallowed. • Irritating to eyes and skin. • Do not breathe vapors and fumes • Wear suitable protective clothing and gloves. • In case of contact with eyes immediately wash with plenty of water and seek medical advice. • If swallowed, seek medical advice and contact poison control center.

This product is intended solely for research use only. Full material safety data sheet document is available at www.lunanano.com.

For further information about this and other products please visit our website at www.lunanano.com.