

ECOCALCIUM



Calcium Alginate Dressing

Highly Absorbent, Maintain Moist Environment

ISO
9001:2015
CERTIFIED

Calcium Alginate Dressing

RECOMMENDED USE

- ✓ All wound depths
- ✓ Moderate/heavy drainage
- ✓ As a primary dressing

RECOMMENDED SECONDARY DRESSINGS

- ✓ Foam dressing
- ✓ Bordered gauze dressing

INDICATIONS

- ✓ Pressure injuries
- ✓ Partial & full-thickness wounds
- ✓ Leg ulcers
- ✓ Diabetic ulcers
- ✓ Surgical wounds
- ✓ Donor sites
- ✓ Lacerations and abrasions
- ✓ First- and second-degree burns

CHANGE FREQUENCY

- ✓ May be left in place for up to 7 days
- ✓ Dressing change frequency will depend on the amount of drainage

CONTRAINDICATIONS

- ✓ Third-degree burns
- ✓ To control heavy bleeding
- ✓ As a surgical sponge
- ✓ Dry or lightly draining wounds
- ✓ Individuals with a known sensitivity to alginates



About EcoCalcium

- ✓ Sterile
- ✓ Single use only
- ✓ Easy to remove
- ✓ Highly absorbent
- ✓ Minimal fibre shed
- ✓ Moist wound environment
- ✓ Conforms to wound contours

Product Code	Size	Packing	GTIN
181-01	5 x 5cm	20boxes/ctn	9337244008440
181-02	10 x 10cm	20boxes/ctn	9337244008457

