

Enclosed Filtration System
BECO INTEGRA® LAB 220 SP



Sheet Filtration System for Single and Two Stage Filtration

BECO INTEGRA LAB 220 SP from Eaton's Begerow Product Line, consisting of a feed pump and filter, was developed for processing small volumes in the laboratory, and for small productions batches an upscaling in the food, chemical, and cosmetic industries

Specific advantages of BECO INTEGRA LAB 220 SP:

- Enclosed filtration system
- Single and two-stage filtration
- Portable unit with centrifugal pump

BECO INTEGRA LAB 220 IP utilizes BECO® depth filter sheets for a nominal pore size from 17 – 0.1 µm.



Technical Data

Max. permissible operating pressure:	58 psi (400 kPa/4 bar)
Max. permissible differential pressure:	43.5 psi (300 kPa/3 bar)
Max. permissible operating temperature:	176 °F (80 °C)

Filter

Filter sheet format:	8.74" (222 mm) round
Number of filter sheets	
Single stage filtration:	2 sheets
Two stage filtrations:	1 x 2 or 2 x 2 sheets
Active filter area:	0.33 ft² (0.03 m²)/ filter sheet
Cake frame width:	0.39" (10 mm)
Max. cake volume:	0.08 gal (0.3 l)

Housing Material

Product contacting parts:	Stainless Steel AISI 304
Gaskets:	NBR



Feed Pump

Parts contacting liquid/product:	Stainless steel
Sealing of housing:	NBR
Rotation speed:	2,800 rpm
Suction height:	Max. 59" (1.5 m) SL
Frequency:	50 Hz
Voltage:	220/230 Volt
Power:	0.37 kW
Protection class:	IP 44

Cable with on-off switch, approx. 78.7" (2 m)

Description	Article No.
BECO INTEGRA LAB 220 SP (centrifugal pump)	59262104
Spare control leakage shaft seal for centrifugal pump	59262110
O-ring sealing pump housing, NBR	P7800314
Set of hoses, incl. connecting pieces in SS	P7800110
Lever for round filter	P7800105
Star knob 50 x M10 (knurled nut)	P7800101

Description	EPDM	Silicone	Viton	FEP coated
O-ring round filter BECO INTEGRA LAB 220	59262117	P7800302	59262118	P7800304



North America - HQ

44 Apple Street,
Tinton Falls, NJ 07724
Toll Free: 800 656-3344
(North America Only)
Voice: +1 732 212-4700

Europe/Africa/Middle East

Auf der Heide 2
53947 Nettersheim, Germany
Voice: +49 2486 809-0

Internormen Product Line

Friedensstraße 41
68804 Altlussheim, Germany
Voice: +49 6205 2094-0

Begerow Product Line

An den Nahewiesen 24
55450 Langenlonsheim, Germany
Voice: +49 6704 204-0

Brazil

Av. Julia Gaioli, 474 - Bonsucesso
07251-500 - Guarulhos
Brazil
Voice: +55 11 2465 8822

China

No. 7 Lane 280 Linhong Road,
Changning District,
Shanghai 200335, China
Voice: +86 21 5200 0422

Singapore

4 Loyang Lane #04-01/02
Singapore 508914
Voice: +65 6825 1668

For more information, please e-mail us at filtration@eaton.com or visit us online at eaton.com/filtration for a complete list of Eaton's filtration products.

Not all products in Eaton's Begerow Product Line are available in all regions. Please contact your local Eaton office to determine availability.

© 2012 Eaton Corporation. All Rights Reserved. All trademarks and registered trademarks are the property of their respective owners.

All information and recommendations appearing in this brochure concerning the use of products described herein are based on tests believed to be reliable. However, it is the user's responsibility to determine the suitability for his own use of such products. Since the actual use by others is beyond our control, no guarantee, expressed or implied, is made by Eaton as to the effects of such use or the results to be obtained. Eaton assumes no liability arising out of the use by others of such products. Nor is the information herein to be construed as absolutely complete, since additional information may be necessary or desirable when particular or exceptional conditions or circumstances exist or because of applicable laws or government regulations.