

ALL STAR FOUNT 2000 TE

Fountain Solution

Meets FOGRA Specs for Sheet-fed Applications

(Heidelberg –Man Roland- Koenig & Bauer)

All Star Fount 2000 TE while based on our most industry proven fountain solution technology, it contains much higher amounts of film formers, lubricants which help keep plates wetter, blankets cleaner and roll ups fast. Buffered to keep its strength to help assure quality is maintained.

- Runs with all offset inks (Conventional, UV, HUV, Hybrid).
- Excellent wetting and viscosity building properties, keeps dynamic surface tension low (lows 30s).
- Helps keep reverses razor sharp, cleaner screens, improved ink trapping and denser solids.
- Compatible with all metal plates, conventional, CTP and non-metal laser plates.

Direction For Use:

Recommended dosage range - 3 oz. to 5 oz. per gallon.
Conductivity per ounce is 200 micromhos.
We recommend starting with 4.5 oz./gallon.
Once running, adjust dosage up or down as needed.

Specifications:

VOC content lb. / gallon	3.3
VOC content grams / liter	390
Flashpoint	None
Color	Colorless
CA. Compliant	Yes



No. 31017



Recommended dosage range 2.5% to 4%.
Conductivity per 1% is 256 micromhos.
We recommend starting with 3%, adjusting dosage up or down as needed.

Available In:

Product Code	Unit	Case Qty.
10653D	1 Gallon (3.785 ltr)	4
10653E	5 Gallons (18.9 ltr)	Each
10653G	55 Gallons (208 ltr)	Each
10653T	275 Gallons (1040 ltr)	Each



VOC COMPLIANT Air Quality Management District Compliant

